

# Medtronic

## SpiderFX™

### Embolic Protection Device

### Capture What Matters.

#### Guidewire of Choice

- Works with 0.014-inch or 0.018-inch guidewire to cross challenging carotid lesions.

#### Extensive Portfolio

- Treat vessels from 3 mm to 7 mm with a variety of sizes; available for carotid interventions;

#### Filter Stability

- Controlled filter positioning throughout the intervention and during device exchanges.
- Braided nitinol design provides full-wall apposition.
- Capture wire designed to rotate and move longitudinally independent of the filter.

#### Enhanced Visibility

- Clearly visible radiopaque markers and direct mouth indicator enable quick and controlled positioning of the filter throughout the intervention.

3.0 mm    4.0 mm    5.0 mm    6.0 mm    7.0 mm

